

German Research Products - GRP GmbH

In der Stockwiese 26

D-85410 Haag/Amper, Germany

Email: info@grp-ak.de Phone: +49 (0)8167 6335

Product Datasheet

Goat anti-Llama IgG (H&L), TRITC conjugated GRP12514

Species/Host Goat

Immunogen Purified llama IgG (H&L) AAQ19986

Form/Appearance Lyophilized

Storage Store lyophilized material at 2-8°C.For long time storage after

reconstitution, dilute the antibody solution with glycerol to a final concentration of 50% glycerol and store as liquid at -20°C, to prevent loss of enzymatic activity. For example, if you have reconstituted 1 mg of antibody in 1.1 ml of sterile water add 1.1 ml of glycerol. Such solution will not freeze in -20°C. If you are using a 1:5000 dilution prior to diluting with glycerol, then you would need to use a 1:2500 dilution after adding glycerol. Prepare working dilution prior to use and then discard. Be sure to mix well

but without foaming.

Note For research use only.

Clonality Polyclonal

Purity Affinity purified IgG

Dilution Range The optimal working dilution should be determined by the

investigator. Suggested starting dilution(s): 1:20-1:2000 for

most applications

Application Notes Additional Information: Purity of this preparation is > 95% based on SDS-PAGE.

Antibody concentration is 1.0 mg/ml (E 1% at 280 nm = 13.0). Antibody is supplied in 10 mM sodium phosphate, 0.15 M sodium chloride, pH 7.2.1 % (w/v) B, Protease/lgG free. Contains 0.05% (w/v) sodium azide as preservative of bacterial growth.Based on immunoelectrophoresis, this antibody reacts with:heavy chains on llama lgG, light chains on all llama immunoglobulins.

Based on immunoelectrophoresis, no reactivity is observed to:

non-immunoglobulin llama serum proteins. This antibody will react with VHH of llama IgG's. Background: Goat anti-llama IgG (H&L) is a secondary antibody

conjugated to TRITC, which binds to llama immunoglobulins (H&L) in

immunological assays. Antibody is affinity purified using solid phase llama IgG (H&L). Tetramethylrhodamine-5-isothiocyanate (TRITC) has Amax = 550 nm; Emax = 570 nm. Reconstitution: For reconstitution add 1.1 ml of sterile water. Let it stand 30 minutes at room temperature to dissolve. Prepare fresh working

dilutions daily.