


Product Datasheet

SWINE IgG (H&L) (RABBIT) Antibody Peroxidase Conjugated GRP5289

Description	Anti-SWINE IgG (H&L) (RABBIT) Antibody Peroxidase Conjugated
Species/Host	Rabbit
Reactivity	Pig
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Form/Appearance	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	2.0 mg/mL
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Source	Anti-Swine IgG whole molecule was produced by repeated immunization with Swine IgG whole molecule in rabbit.
Dilution Range	ELISA: 1:50,000 - 1:100,000, Immunochemistry: 1:500 - 1:2,500, Western Blot: 1:1,000 - 1:5,000
Application Notes	Suitable for immunoblotting (western or dot blot), ELISA, immunoperoxidase electron microscopy and immunohistochemistry as well as other peroxidase-antibody based enzymatic assays requiring lot-to-lot consistency.



ELISA Results of Rabbit Anti-Swine IgG Antibody Peroxidase Conjugate