


Product Datasheet

BIOTIN PEROXIDASE Conjugated GRP10327

Description	BIOTIN PEROXIDASE Conjugated
Conjugation	Biotin
Tested Applications	ELISA, IHC-P, WB
Form/Appearance	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	1.0 mg/mL
Note	For research use only.
Purity	Biotin Peroxidase Conjugated was prepared from chromatographically purified biotin. Biotin Peroxidase Conjugated was assayed by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase and anti-Biotin.
Dilution Range	ELISA: 1:20,000 - 1:200,000, Immunochemistry: 1:1,000 - 1:5,000, Western Blot: 1:10,000 - 1:40,000
Application Notes	Biotin Peroxidase Conjugated can be utilized in ELISA and Western Blotting experiments where the assay's target of interest is coupled with streptavidin.



Biotin Peroxidase
Conjugated -
Western Blot