

German Research Products - GRP GmbH

In der Stockwiese 26

D-85410 Haag/Amper, Germany

Email: info@grp-ak.de Phone: +49 (0)8167 6335

Product Datasheet

Goat anti-Rabbit IgG (H&L), DyLight® 550 conjugated GRP12776

Species/Host Goat

Reactivity Rabbit

Predicted Reactivity Rabbit IgG Heavy and Light chains (H&L)

Tested Applications ICC, IF, IHC

Immunogen Purified Rabbit IgG, whole molecule

Form/Appearance Lyophilized

Storage Store lyophilized material at 2-8°C. Product is stable for 4 weeks

at 2-8°C after rehydration. For long time storage after

reconstitution, dilute the antibody solution with glycerol to a final concentration of 50% glycerol and store as liquid at -20°C, to prevent loss of enzymatic activity. For example, if you have reconstituted 1 mg of antibody in 1.1 ml of sterile water add 1.1 ml of glycerol. Such solution will not freeze in -20°C. If you are using a 1:5000 dilution prior to diluting with glycerol, then you would need to use a 1:2500 dilution after adding glycerol. Prepare working dilution prior to use and then discard. Be sure to mix well

but without foaming.

Note For research use only.

Clonality Polyclonal

Purity Affinity purified goat IgG

Dilution Range 1 : 50- 1 : 5 000 (ICC), 1 : 20- 1 : 2000 (IHC), 1 : 3000 (IF)

Application Notes Additional Information: Based in immunoelectrophoresis, this antibody reacts

with heavy chains on rabbit IgG and light chains on all rabbit

immunoglobulins.No reactivity is observed to non-immunoglobulin rabbit serum proteins based in immunoelectrophoresis.Purity of this antibody is > 95% based on SDS-PAGE. Conjugate is present in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/lgG free. 0.05 % (w/v) sodium azide is added as preservative.DyLight® 550 (Ex = 550 nm; Em = 576 nm) Background:

Goat anti-rabbit IgG, DyLight® 550 Conjugate is a secondary antibody conjugated toDyLight® 550, which binds to all rabbit IgGs in immunological assays.DyLight® is a trademark of Thermo Fisher Scientific, Inc. and its

subsidiaries. Reconstitution: For reconstitution add 1.1 ml of sterile water. Let it

stand 30 minutes at room temperature to dissolve. Prepare fresh working

dilutions daily.