

Product Datasheet

HORSE IgG F(ab')₂ (RABBIT) Antibody Peroxidase Conjugated GRP4470

Description Anti-HORSE IgG F(ab')₂ (RABBIT) Antibody Peroxidase Conjugated

Species/Host Rabbit

Conjugation HRP

Tested Applications ELISA, IHC-P, WB

Form/Appearance 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Concentration 1.5 mg/mL

Note For research use only.


Isotype IgG

Clonality Polyclonal

Source Horse IgG F(ab')₂ fragment

Dilution Range ELISA: 1:10,000 - 1:50,000, Immunochemistry: 1:500 - 1:2,500, Western Blot: 1:1,000 - 1:5,000

Application Notes Anti-Horse IgG F(ab')₂ Peroxidase conjugate is suitable for immunoblotting (western or dot blot), ELISA, immunoelectron microscopy and immunohistochemistry as well as other antibody-based enzymatic assays requiring lot-to-lot consistency.



ELISA Results of Rabbit Anti-Horse IgG F(ab')₂ Antibody Peroxidase Conjugated